

KN95 MASK RESELLER ACKNOWLEDGMENT

I acknowledge that I have been provided a document with the Product Information, Instructions and Warnings below, and I agree to provide a copy of such document to any person, business or facility of any kind to which I provide or cause to be provided this Face Mask. I further agree to advise such person, business or facility to prominently display a copy of these Instructions and Warnings in locations that are visible and accessible to users of these Face Masks, and to otherwise provide users with a copy of the Product Information, Instructions and Warnings.

Name and Title *(this serves as electronic signature)* Company Date

PRODUCT INFORMATION, INSTRUCTIONS FOR USE, AND WARNINGS

Mask Model No.: KN95

Manufactured by: Shenzhen Zhiyou Medical Equipment Co., Ltd.
3rd floor, 4th building, Transcend Phase II Industrial Zone, Yingrenshi
Community, Shiyan Street, Baoan District,
Shenzhen, Guangdong, 518108, China

Distributed by: The Millcraft Paper Company
6800 Grant Avenue
Cleveland, Ohio 44105
www.millcraft.com *(for terms and conditions)*

Intended Use: This KN95 mask is a non-sterile, single use, non-powered particle-filtering half mask that covers the nose, mouth and chin, and when properly fit and worn, is intended to provide adequate sealing on the face of the wearer against particulate contaminants in the ambient air. This KN95 mask has particulate filtering efficiency (PFE) of $\geq 95\%$ for particles 0.3 microns in diameter or larger.

 **WARNINGS**

- **Achieving an adequate seal to the face when using this respirator mask (“mask) is essential. User should conduct a seal check prior to use. Facial hair such as beards, sideburns, or some mustaches, may interfere with the proper fit and performance of these masks.**
- **These masks are NOT FDA cleared or approved but HAVE been given Emergency Use Authorization by the FDA because of shortages resulting from the Coronavirus Disease 2019 (COVID-19) outbreak.**
- **These masks have NOT been determined to be safe or effective in the prevention of COVID-19 transmission. It is the sole responsibility of the user to determine the proper and best use for these masks.**
- **We make no claims regarding effectiveness; antimicrobial or antiviral protection; infection prevention or reduction; flammability; or any other claims with respect to the performance of these masks.**
- **These masks are non-sterile and are NOT surgical masks. They are not suitable for use in surgical or other sterile environments or where significant exposure to bodily liquids or other hazardous fluids may be expected.**
- **These masks are not suitable for use in a clinical setting where the infection risk level through inhalation exposure is high. These masks are not suitable for use in the presence of a high intensity heat source or flammable gas.**
- **This mask is made of non-woven fabric and melt-blown cloth (80% polypropylene, 20% ES fiber).**