

3 PLY FACE MASK RESELLER ACKNOWLEDGMENT



I acknowledge that I have been provided a document with the Product Information, Instructions and Warnings below, and I agree to provide a copy of such document to any person, business or facility of any kind to which I provide or cause to be provided this Face Mask. I further agree to advise such person, business or facility to prominently display a copy of these Instructions and Warnings in locations that are visible and accessible to users of these Face Masks, and to otherwise provide users with a copy of the Product Information, Instructions and Warnings.

Name and Title *(this serves as electronic signature)*

Company

Date

PRODUCT INFORMATION, INSTRUCTIONS, AND WARNINGS

Product: 3 Ply Face Mask

Manufactured by: Shandong Bossen Medicine Engineering Technology Co LTD Floor 2, Block B, Qilu Science and Technology Incubators District Qihe Country Dezhou, Shandong, 253000 China

Distributed by: The Millcraft Paper Company
6800 Grant Avenue
Cleveland, Ohio 44105
www.millcraft.com *(for terms and conditions)*

Intended Use: This face mask is a non-sterile, single use, mask intended for use by members of the general public, including healthcare workers in healthcare settings, as personal protective equipment, during the duration of the COVID-19 pandemic.



WARNINGS

- This face mask is **NOT** a surgical mask and has **NOT** been FDA cleared or approved.
- This face mask **HAS** been authorized by FDA under an Emergency Use Authorization (EUA) for use by healthcare providers as personal protective equipment to help prevent the spread of infection or illness in healthcare settings and by the general public to help slow the spread of the virus during the COVID-19 pandemic.
- This face mask is latex free and made of non-woven, 3-layer polypropylene and latex free elastic
- This face mask is **NOT** recommended for use in any surgical setting or where significant exposure to liquid, bodily or other hazardous fluids may be expected.
- This face mask is **NOT** recommended for in a clinical setting where the infection risk level through inhalation exposure is high.
- **Do NOT** use this mask for antimicrobial or antiviral protection, for prevention or reduction of infection or for particulate filtration.
- This face mask has **NOT** been determined to be safe or effective in the prevention of COVID-19 transmission. It is the sole responsibility of the user to determine the proper and best use for these masks.
- **Do NOT** use this mask in the presence of a high intensity heat source or flammable gas.